



**Task
48**

**Determination of the Planck constant using
electroluminescent diodes**

Instruction

1. The list of necessary equipment

- Tunable power supply
- Electroluminescent diode
- Digital multimeters
- Monochromator
- Photoresistor

2. Goals

- Measurement of the current-voltage characteristic for different diodes in forward direction (forward bias)
- Determination of the emission wavelength for investigated diodes
- Calculation of Planck constant and related uncertainties.

3. Experimental setup




Fig. 1. The electrical circuit used to measure the current-voltage characteristic for forward-biased diodes.




Fig. 2. The block-diagram of the experimental setup used for determination of the current-voltage characteristic of a LED diode.

4. Measurements

- a) Connect the required elements according to Fig. 2.
- b) Set the power supply potentiometer to the lowest position (turn potentiometer to the left).
- c) Set the power supply to forward bias.
- d) Measure the current-voltage characteristic. Use the potentiometer to increase the applied voltage. Conduct these measurements in as wide range of voltages as possible. Write down the measured voltage U and current I . You should have at least 20 data points measured.
- e) Determine the emission wavelength of a diode using one of the following methods:
 - i. By eye:
 - Set the power supply to maximum.
 - Place the diode at the rear entrance of the monochromator.
 - Change the wavelength of light transmitted to the exit slit of the monochromator (by turning the selection knob) so that its intensity is the highest.
 - Read the wavelength λ from the selection knob.
 - Assume $\Delta\lambda = 10 \text{ nm}$ for the uncertainty.
 - ii. Photoresistor:
 - Set the power supply to maximum.
 - Place the diode at the rear entrance of the monochromator.
 - Place the photoresistor in front of the monochromator.
 - Change the wavelength of light transmitted to the exit slit of the monochromator (by turning the selection knob) so that the resistance magnitude is the lowest.
 - Read the wavelength λ from the selection knob.
 - Assume $\Delta\lambda = 5 \text{ nm}$ for the uncertainty.
- f) Repeat the measurements for other diodes.

5. Calculations for report

- a) Plot current-voltage characteristics for measured diodes.
- b) Using linear regression method, fit a straight line $y = ax + b$ to the linear part of the measured current-voltage characteristic (forward direction only).
- c) Calculate the potential barrier U_b of the diode, according to equation: $U_b = -\frac{b}{a}$.
- d) Estimate the uncertainty ΔU_b related to U_b
- e) Calculate the Planck constant using the following equation

$$h = \frac{e}{c} \lambda U_b$$

where e is the elementary charge and c is the velocity of light in vacuum.

- f) Calculate the uncertainty Δh related to the determined Planck constant.

6. Table for measurements of the current-voltage characteristic

Table 1. Results of your experiment.

lp.	U [V]	$u(U)$ [V]	I [mA]	$u(I)$ [mA]	λ [nm]	$u(\lambda)$ [nm]	U_b [V]	$u_c(U_b)$ [V]	h [J·s]	$u_c(h)$ [J·s]
1										
2										
3										
4										
5										
6										
7										
8										
9										
10										
11										
12										
13										
14										
15										
16										
17										
18										
19										
20										